

## Human DNASE1 (DNASE1) Sandwich ELISA Kit (STJE0000631)

STJE0000631

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Kits   |
| <b>Short Description</b> | This DNASE1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA  |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                         |   |
|-------------------------|---|
| <b>Detection</b>        | Colorimetric  |
| <b>Sensitivity</b>      | 0.258 ng/mL   |
| <b>Detection Limit</b>  | 0.78-50 ng/mL   |
| <b>Assay Time</b>       | 4.0 hours   |
| <b>Storage</b>          | Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert. |
| <b>Instruction Note</b> | FOR RESEARCH USE ONLY (RUO).  |

### TARGET INFORMATION

|                     |             |
|---------------------|-------------|
| <b>Sample Types</b> |             |
| <b>Gene</b>         | DNASE1      |
| <b>Symbol</b>       |             |
| <b>Uniprot ID</b>   | DNAS1_HUMAN |