

Human FNDC5 (FNDC5) Sandwich ELISA Kit (STJE0000767)

STJE0000767

GENERAL INFORMATION

Product Type	Kits
Short Description	This FNDC5 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	14.51pg/mL
Detection Limit	2000-31.2pg/ml
Assay Time	
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	FNDC5
Symbol	
Uniprot ID	FNDC5_HUMAN