

Human HIF1a (HIF1a) Sandwich ELISA Kit (STJE0000920)

STJE0000920

GENERAL INFORMATION

Product Type	Kits
Short Description	This HIF1a Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	35.1pg/ml
Detection Limit	78.1-5000pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Gene	HIF1A
Symbol	
Uniprot ID	HIF1A_HUMAN