

Human ENO2 (ENO2) Sandwich ELISA Kit (STJE0001313)

STJE0001313

GENERAL INFORMATION

Product Type	Kits
Short Description	This ENO2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.17 ng/mL
Detection Limit	20-0.312ng/ml
Assay Time	3h
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	ENO2
Symbol	
Uniprot ID	ENOG_HUMAN