

## Human PD-L1/B7-H1/CD274 (PD-L1/B7-H1/CD274) Sandwich ELISA Kit (STJE0001891)

STJE0001891

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This PD-L1/B7-H1/CD274 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Colorimetric
<b>Sensitivity</b>	16.29 pg/mL
<b>Detection Limit</b>	156-10000pg/ml
<b>Assay Time</b>	
<b>Storage</b>	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
<b>Instruction</b>	
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

**Sample Types**  
**Gene**  
**Symbol**  
**Uniprot ID**