

Human TNFSF9 (TNFSF9) Sandwich ELISA Kit (STJE0001893)

STJE0001893

GENERAL INFORMATION

Product Type	Kits
Short Description	This TNFSF9 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	10.25pg/mL
Detection Limit	2000-31.2pg/ml
Assay Time	
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	TNFSF9
Symbol	
Uniprot ID	TNFL9_HUMAN