

Human DKK3 (DKK3) Sandwich ELISA Kit (STJE0001898)

STJE0001898

GENERAL INFORMATION

Product Type	Kits
Short Description	This DKK3 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	45.3 pg/ml
Detection Limit	125-8000 pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, cell culture supernates
Gene	DKK3
Symbol	
Uniprot ID	DKK3_HUMAN