

## SARS-CoV-2 S1+S2 ECD Sandwich ELISA Kit (STJE0001904)

STJE0001904

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This SARS-CoV-2 S1+S2 ECD (S-ECD) Protein Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in sars-cov-2 cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	SARS-CoV-2

### PRODUCT PROPERTIES

<b>Detection</b>	Colorimetric
<b>Sensitivity</b>	
<b>Detection Limit</b>	
<b>Assay Time</b>	
<b>Storage</b>	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
<b>Instruction</b>	
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

**Sample Types**  
**Gene Symbol**  
**Uniprot ID**