

## Bovine, Horse, Sheep E2 (Estradiol) Competitive ELISA Kit (STJE0001908)

STJE0001908

### GENERAL INFORMATION


|                          |                                                                                                                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                                                       |
| <b>Short Description</b> | This E2 Competitive ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in bovine, horse, sheep cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA                                                                                                                                                                                                      |
| <b>Reactivity</b>        | Bovine, Horse, Sheep                                                                                                                                                                                       |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                                             |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | Colormetric                                                                                                                                                                                 |
| <b>Sensitivity</b>     | 5.63pg/mL                                                                                                                                                                                   |
| <b>Detection Limit</b> | 9.37~600pg/mL                                                                                                                                                                               |
| <b>Assay Time</b>      | 2.5h                                                                                                                                                                                        |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Instruction</b>     |                                                                                                                                                                                             |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).                                                                                                                                                                |

### TARGET INFORMATION

|                     |                  |
|---------------------|------------------|
| <b>Sample Types</b> | Serum and plasma |
| <b>Gene</b>         |                  |
| <b>Symbol</b>       |                  |
| <b>Uniprot ID</b>   |                  |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081