

Human IGFBP4 (Insulin-like growth factor-binding protein 4) Sandwich ELISA Kit (STJE0002234)

STJE0002234

GENERAL INFORMATION

Product Type	Kits
Short Description	This IGFBP4 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	1.88ng/mL
Detection Limit	3.13-200ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	IGFBP4
Uniprot ID	IBP4_HUMAN