

Human TNFRSF13C (Tumor necrosis factor receptor superfamily member 13C) Sandwich ELISA Kit (STJE0002307)

STJE0002307

GENERAL INFORMATION

Product Type	Kits
Short Description	This TNFRSF13C Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	0.09ng/mL
Detection Limit	0.16~10ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	TNFRSF13C
Symbol	
Uniprot ID	TR13C_HUMAN