

Human GPS1 (COP9 signalosome complex subunit 1) Sandwich ELISA Kit (STJE0003741)

STJE0003741

GENERAL INFORMATION

Product Type	Kits
Short Description	This GPS1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	46.88ng/mL
Detection Limit	78.13~5000ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	GPS1
Symbol	
Uniprot ID	CSN1_HUMAN