

Human CXCL16 (C-X-C motif chemokine 16) Sandwich ELISA Kit (STJE0004021)

STJE0004021

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This CXCL16 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|---|
| Detection | Colorimetric |
| Sensitivity | 0.19ng/mL |
| Detection Limit | 0.31~20ng/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene | CXCL16 |
| Symbol | |
| Uniprot ID | CXL16_HUMAN |