

Human CEACAM5 (Carcinoembryonic antigen-related cell adhesion molecule 5) Sandwich ELISA Kit (STJE0004059)

STJE0004059

GENERAL INFORMATION

Product Type	Kits
Short Description	This CEACAM5 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.1ng/mL
Detection Limit	0.16-10ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	CEACAM5
Symbol	
Uniprot ID	CEAM5_HUMAN