

Human C4B; C4B_2 (Complement C4-B) Sandwich ELISA Kit (STJE0004071)

STJE0004071

GENERAL INFORMATION

Product Type	Kits
Short Description	This C4B; C4B_2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.94ng/mL
Detection Limit	1.56~100ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	C4B.C4B_2
Symbol	
Uniprot ID	CO4B_HUMAN