

Human ORM2 (Alpha-1-acid glycoprotein 2) Sandwich ELISA Kit (STJE0004098)

STJE0004098

GENERAL INFORMATION

Product Type	Kits
Short Description	This ORM2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	0.19ng/mL
Detection Limit	0.31~20ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	ORM2
Symbol	
Uniprot ID	A1AG2_HUMAN