

Human CDH5 (Cadherin-5) Sandwich ELISA Kit (STJE0004115)

STJE0004115

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This CDH5 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|---|
| Detection | Colorimetric |
| Sensitivity | 18.75pg/mL |
| Detection Limit | 31.25~2000pg/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene | CDH5 |
| Symbol | |
| Uniprot ID | CADH5_HUMAN |