

Mouse CD257/Tnfsf13b (Tumor necrosis factor ligand superfamily member 13B) Sandwich ELISA Kit (STJE0004205)

STJE0004205

GENERAL INFORMATION

Product Type	Kits
Short Description	This CD257/Tnfsf13b Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	37.5pg/mL
Detection Limit	62.5~4000pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	Tnfsf13b
Symbol	
Uniprot ID	TN13B_MOUSE