

Mouse Rarres2 (Retinoic acid receptor responder protein 2) Sandwich ELISA Kit (STJE0004212)

STJE0004212

GENERAL INFORMATION

Product Type	Kits
Short Description	This Rarres2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	28.13pg/mL
Detection Limit	46.88~3000pg/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	Rarres2
Uniprot ID	RARR2_MOUSE