

Mouse Cxcl13 (C-X-C motif chemokine 13) Sandwich ELISA Kit (STJE0004328)

STJE0004328

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This Cxcl13 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|---|
| Detection | Colometric |
| Sensitivity | 9.38pg/mL |
| Detection Limit | 15.63~1000pg/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene | Cxcl13 |
| Symbol | |
| Uniprot ID | CXL13_MOUSE |