

Mouse Ggt1 (Glutathione hydrolase 1 proenzyme) Sandwich ELISA Kit (STJE0004590)

STJE0004590

GENERAL INFORMATION


| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Kits |
| Short Description | This Ggt1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detection | Colometric |
| Sensitivity | 0.09ng/mL |
| Detection Limit | 0.16~10ng/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|-------------------------------------------|
| Sample Types | Serum, plasma and other biological fluids |
| Gene | Ggt1 |
| Symbol | |
| Uniprot ID | GGT1_MOUSE |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081