

Mouse Fabp3 (Fatty acid-binding protein, heart) Sandwich ELISA Kit (STJE0005169)

STJE0005169

GENERAL INFORMATION

Product Type	Kits
Short Description	This Fabp3 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	37.5pg/mL
Detection Limit	62.5~4000pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	Fabp3
Symbol	
Uniprot ID	FABPH_MOUSE