

Mouse Igfbp5 (Insulin-like growth factor-binding protein 5) Sandwich ELISA Kit (STJE0005224)

STJE0005224

GENERAL INFORMATION

Product Type	Kits
Short Description	This Igfbp5 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	75pg/mL
Detection Limit	125~8000pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	Igfbp5
Symbol	
Uniprot ID	IBP5_MOUSE