

## Mouse Ighm (Immunoglobulin heavy constant mu) Sandwich ELISA Kit (STJE0005352)

STJE0005352

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Kits   |
| <b>Short Description</b> | This Ighm Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA  |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                        |   |
|------------------------|---|
| <b>Detection</b>       | Colometric  |
| <b>Sensitivity</b>     | 1.88ng/mL   |
| <b>Detection Limit</b> | 3.13~200ng/mL   |
| <b>Assay Time</b>      | 3.5h  |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Instruction</b>     |   |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).  |

### TARGET INFORMATION

|                     |   |
|---------------------|---|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene</b>         | Ighm                                      |
| <b>Symbol</b>       |   |
| <b>Uniprot ID</b>   | IGHM_MOUSE                                |