

Mouse Cxcl12 (Stromal cell-derived factor 1) Sandwich ELISA Kit (STJE0005362)

STJE0005362

GENERAL INFORMATION

Product Type	Kits
Short Description	This Cxcl12 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	46.88pg/mL
Detection Limit	78.13~5000pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	Cxcl12
Symbol	
Uniprot ID	SDF1_MOUSE