

Rat Col1a1 (Collagen alpha-1 (I) chain) Sandwich ELISA Kit (STJE0005452)

STJE0005452

GENERAL INFORMATION

Product Type	Kits
Short Description	This Col1a1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.19ng/mL
Detection Limit	0.31~20ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	Col1a1
Symbol	
Uniprot ID	CO1A1_RAT