

## Rat Cp (Ceruloplasmin) Sandwich ELISA Kit (STJE0005558)

STJE0005558

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Kits   |
| <b>Short Description</b> | This Cp Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA  |
| <b>Reactivity</b>        | Rat  |

### PRODUCT PROPERTIES

|                        |   |
|------------------------|---|
| <b>Detection</b>       | Colormetric   |
| <b>Sensitivity</b>     | 3.75ng/mL   |
| <b>Detection Limit</b> | 6.25~400ng/mL   |
| <b>Assay Time</b>      | 3.5h  |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Instruction</b>     |   |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).  |

### TARGET INFORMATION

|                     |   |
|---------------------|---|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  | Cp  |
| <b>Uniprot ID</b>   | CERU_RAT                                  |