

Rat Tf (Serotransferrin) Sandwich ELISA Kit (STJE0006563)

STJE0006563

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This Tf Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|---|
| Detection | Colormetric |
| Sensitivity | 37.5pg/mL |
| Detection Limit | 62.5~4000pg/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene Symbol | Tf |
| Uniprot ID | TRFE_RAT |