

BH4 (Tetrahydrobiopterin) Competitive ELISA Kit (STJE0006619)

STJE0006619

GENERAL INFORMATION


Product Type	Kits
Short Description	This BH4 Competitive ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in all cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	All

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	18.75pg/mL
Detection Limit	31.25~2000pg/mL
Assay Time	
Storage	
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Cell supernatant
Gene	
Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081