

Sars-Cov-2 SARS-CoV-2 (SARS-CoV-2 Spike Protein IgG) Indirect ELISA Kit (STJE0006696)

STJE0006696

GENERAL INFORMATION

Product Type	Kits
Short Description	This SARS-CoV-2 Indirect ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in sars-cov-2 cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	
Detection Limit	
Assay Time	1.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum and plasma
Gene	
Symbol	
Uniprot ID	