

## Sars-Cov-2 SARS-CoV-2 (SARS-CoV-2 Spike Protein IgM) Indirect ELISA Kit (STJE0006697)

STJE0006697

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This SARS-CoV-2 Indirect ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in sars-cov-2 cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	SARS-CoV-2

### PRODUCT PROPERTIES

<b>Detection</b>	Colometric
<b>Sensitivity</b>	
<b>Detection Limit</b>	
<b>Assay Time</b>	1.5h
<b>Storage</b>	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
<b>Instruction</b>	
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	Serum and plasma
<b>Gene</b>	
<b>Symbol</b>	
<b>Uniprot ID</b>	