

SARS-CoV-2 Spike Protein IgM ELISA Kit (STJE0006697)

STJE0006697

GENERAL INFORMATION

Product Type	Kits
Short Description	This SARS-CoV-2 Indirect ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	
Detection Limit	
Assay Time	
Storage	
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma
Gene Symbol	
Uniprot ID	