

## Sars-Cov-2 SARS-CoV-2 (SARS-CoV-2 Nucleocapsid Protein) Sandwich ELISA Kit (STJE0006698)

STJE0006698

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This SARS-CoV-2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in sars-cov-2 cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	SARS-CoV-2

### PRODUCT PROPERTIES

<b>Detection</b>	Colometric
<b>Sensitivity</b>	0.12ng/mL
<b>Detection Limit</b>	0.39~25ng/mL
<b>Assay Time</b>	3.5h
<b>Storage Instruction</b>	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma and other biological fluids
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	