

## Human SARS-CoV-2 Neutralization Antibody (SARS-CoV-2 Neutralization Antibody) Competitive ELISA Kit (STJE0006700)

STJE0006700

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This SARS-CoV-2 Neutralization Antibody Competitive ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                        |   |
|------------------------|---|
| <b>Detection</b>       | Colormetric   |
| <b>Sensitivity</b>     |   |
| <b>Detection Limit</b> |   |
| <b>Assay Time</b>      | 1.5h  |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual. |
| <b>Instruction</b>     |   |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).  |

### TARGET INFORMATION

|                     |                  |
|---------------------|------------------|
| <b>Sample Types</b> | Serum and plasma |
| <b>Gene</b>         |                  |
| <b>Symbol</b>       |                  |
| <b>Uniprot ID</b>   |                  |