

Human SARS-CoV-2 Neutralization Antibody (SARS-CoV-2 Neutralization Antibody) Competitive ELISA Kit (STJE0006701)

STJE0006701

GENERAL INFORMATION


Product Type	Kits
Short Description	This SARS-CoV-2 Neutralization Antibody Competitive ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	9.38ng/mL
Detection Limit	15.63~500ng/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081