

Human LCN2 (Neutrophil gelatinase-associated lipocalin) Sandwich ELISA Kit (STJE0006705)

STJE0006705

GENERAL INFORMATION

Product Type	Kits
Short Description	This LCN2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	18.75pg/ml
Detection Limit	31.25~2000 pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	LCN2
Symbol	
Uniprot ID	NGAL_HUMAN