

Human IL6 (Interleukin-6) Sandwich ELISA Kit (STJE0006744)

STJE0006744

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL6 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	1.02pg/mL
Detection Limit	3.13-200pg/mL
Assay Time	1.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum and plasma
Gene	IL6
Symbol	
Uniprot ID	IL6_HUMAN