

Human AFP (Alpha-fetoprotein) Sandwich ELISA Kit (STJE0006756)

STJE0006756

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This AFP Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|---|
| Detection | Colormetric |
| Sensitivity | 0.47ng/mL |
| Detection Limit | 0.78-50ng/mL |
| Assay Time | 1.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|------------------|
| Sample Types | Serum and plasma |
| Gene | AFP |
| Symbol | |
| Uniprot ID | FETA_HUMAN |