

Human CEACAM5 (Carcinoembryonic antigen-related cell adhesion molecule 5) Sandwich ELISA Kit (STJE0006757)

STJE0006757

GENERAL INFORMATION

Product Type	Kits
Short Description	This CEACAM5 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.47ng/mL
Detection Limit	0.78-50ng/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum and plasma
Gene	CEACAM5
Symbol	
Uniprot ID	CEAM5_HUMAN