

Mouse E3 (Estriol) Competitive ELISA Kit (STJE0006759)

STJE0006759

GENERAL INFORMATION

Product Type	Kits
Short Description	This E3 Competitive ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	5.83pg/mL
Detection Limit	15.63-1000pg/mL
Assay Time	1.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum and plasma
Gene	
Symbol	
Uniprot ID	