

Mouse IL1RA (Interleukin 1 Receptor Antagonist) Sandwich ELISA Kit Ready-To-Use (STJE0006850)

STJE0006850

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL1RA Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	31.2-2000pg/mL
Sensitivity	11.2pg/mL
Detection Limit	31.2-2000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	Il1rn
Symbol	
Uniprot ID	IL1RA_MOUSE