

Human IGFBP4 (Insulin Like Growth Factor Binding Protein 4) Sandwich ELISA Kit Ready-To-Use (STJE0006857)

STJE0006857

GENERAL INFORMATION

Product Type	Kits
Short Description	This IGFBP4 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	3.12-200ng/mL
Sensitivity	1.13ng/mL
Detection Limit	3.12-200ng/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma or other biological fluids.
Gene	IGFBP4
Symbol	
Uniprot ID	IBP4_HUMAN