

Mouse TNFRSF5 (Tumor Necrosis Factor Receptor Superfamily, Member 5) Sandwich ELISA Kit Ready-To-Use (STJE0006863)

STJE0006863

GENERAL INFORMATION

Product Type	Kits
Short Description	This TNFRSF5 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	15.6-1000pg/mL
Sensitivity	5.4pg/mL
Detection Limit	15.6-1000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma or other biological fluids.
Gene	Cd40
Symbol	
Uniprot ID	TNR5_MOUSE