

Human EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor) Sandwich ELISA Kit Ready-To-Use (STJE0006889)

STJE0006889

GENERAL INFORMATION

Product Type	Kits
Short Description	This EG-VEGF Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	15.6-1000pg/mL
Sensitivity	7.2pg/mL
Detection Limit	15.6-1000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	PROK1
Symbol	
Uniprot ID	PROK1_HUMAN