

Mouse SELP (Selectin, Platelet) Sandwich ELISA Kit Ready-To-Use (STJE0006925)

STJE0006925

GENERAL INFORMATION

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Kits |
| Short Description | This SELP Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detection | 0.156-10ng/mL |
| Sensitivity | 0.068ng/mL |
| Detection Limit | 0.156-10ng/mL |
| Assay Time | 3 hrs |
| Storage | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|------------------------------------------------------------------------------------------------------|
| Sample Types | serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids. |
| Gene | Selp |
| Symbol | |
| Uniprot ID | LYAM3_MOUSE |