

Mouse VCAM1 (Vascular Cell Adhesion Molecule 1) Sandwich ELISA Kit Ready-To-Use (STJE0006941)

STJE0006941

GENERAL INFORMATION

Product Type	Kits
Short Description	This VCAM1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	0.39-25ng/mL
Sensitivity	0.23ng/mL
Detection Limit	0.39-25ng/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	Vcam1
Symbol	
Uniprot ID	VCAM1_MOUSE