

Mouse TREM1 (Triggering Receptor Expressed On Myeloid Cells 1) Sandwich ELISA Kit Ready-To-Use (STJE0006956)

STJE0006956

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This TREM1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Detection | 15.6-1000pg/mL |
| Sensitivity | 5.4pg/mL |
| Detection Limit | 15.6-1000pg/mL |
| Assay Time | 3 hrs |
| Storage Instruction | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|--|
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Gene | Trem1 |
| Symbol | |
| Uniprot ID | TREM1_MOUSE |