

Human Flt3L (FMS Like Tyrosine Kinase 3 Ligand) Sandwich ELISA Kit Ready-To-Use (STJE0006961)

STJE0006961

GENERAL INFORMATION

Product Type	Kits
Short Description	This Flt3L Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	15.6-1000pg/mL
Sensitivity	7.5pg/mL
Detection Limit	15.6-1000pg/mL
Assay Time	3 hrs
Storage Instruction	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	FLT3LG
Symbol	
Uniprot ID	FLT3L_HUMAN