

Human CTACK (Cutaneous T-Cell Attracting Chemokine) Sandwich ELISA Kit Ready-To-Use (STJE0006997)

STJE0006997

GENERAL INFORMATION

Product Type	Kits
Short Description	This CTACK Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Sensitivity	31.2-2000pg/mL
Detection Limit	10.2pg/mL
Assay Time	31.2-2000pg/mL
Storage Instruction	3 hrs
Note	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene Symbol	CCL27
Uniprot ID	CCL27_HUMAN