

Mouse CLU (Clusterin) Sandwich ELISA Kit Ready-To-Use (STJE0007003)

STJE0007003

GENERAL INFORMATION

Product Type	Kits
Short Description	This CLU Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	18.7-1200ng/mL
Sensitivity	6.1ng/mL
Detection Limit	18.7-1200ng/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, urine, breast milk, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates or other biological fluids.
Gene	Clu
Symbol	
Uniprot ID	CLUS_MOUSE